



Introducing XN and XN-L series:

What is in it  
for you?

Dr Marion Münster  
Manager Medical & scientific Affairs




XN

# XN Concept



XN

Each axis with enhanced value  
– “Modular and Scalable”



# The XN-Series and XN-L Series



XN-350

XN-450

XN-550



XN-1000


XN-2000

XN-9000

XN-Series



# Fluorescence Flow Cytometry





# Low Aspiration Volume



# Reflex Testing






Manual and Automated



# Secure Sample Identification



# Quality Control




# SNCS - Sysmex Network Communication System



- SNCS IQAS Online



## Internal QC



## External QC



- Daily QC measurement of XN Check and XN Check BF



# Reagents: Plug and Play





# Reagent Management



# The XN portrait



Technologies: Fluorescence Flow Cytometry in all modes

Aspiration volume: 88  $\mu$ l in all modes

Throughput: starting from 100 samples/h

Quality control: XN check and XN check BF

Parameters: 28 diagnostic parameters always standard  
*XN-CBC = always with NRBC*  
*XN-DIFF = efficiency of XE-5000*  
16 diagnostic parameters are optional

Configurations: XN-1000, XN-2000, XN-3000 or XN-9000



# XN Concept – Productivity Values

---

**XN**


The letters "XN" in a bold, grey, sans-serif font, positioned above a reflection of the same letters. The background features a series of horizontal lines that create a sense of depth and perspective, suggesting a modern architectural or industrial setting.

## Productivity Values




- *From bench top to automated solution*
- *Suitable for every lab organisation, lab size and workload*
- *Protects your initial investment*
- *Prepared for future growth*

# XN Concept – Productivity Values







# A system for seamless integration of productivity values (= hardware)



A system for seamless integration of clinical values and productivity values coupled with intelligence



# XN Concept – Clinical Values



# XN Concept – Clinical Values

---

**XN**

Adding value...

...to routine haematology testing



**The XN-DIFF:**

**World-leading diagnostic performance and flagging efficiency in all XN modules for true peace-of-mind**

- Rerun/reflex “Low WBC mode” for critically low cell counts
- Multi-dimensional flagging information – highly sensitive identification of pathological cells
- Automated IG count





## Special „Low WBC mode“



- Extended count volume (3 x) – prolonged count, same aspiration volume
- For high accuracy of WBC# and differential in leucopenic samples
- WBC# triggers reflex testing in the “Low WBC mode”  $<1000$  WBC/ $\mu$ L or vote out for Neut#
- Sample is automatically reflexed and rerun in low WBC mode

**Accurate neutrophil count  
when it matters most!**

**“Decide whether to give or  
withhold chemotherapy  
with confidence”**






# Differential count with multi-dimensional abnormal blood cell identification and flagging



## Providing more than just numeric information:

- By special shape recognition, analysis of each subpopulation
- Flexible gating for counting
- Sensitivity comparable with smear morphology



**The presence of abnormal cells will not be missed!**



Pathological cells will not be missed!

- With highly sensitive flagging
- Samples needing manual review are automatically identified
- Technologist is guided in what to look for **BUT**.....

To ensure that microscopic review will provide a report that can be trusted for clinical judgement.....

**Quality of smear and stain MUST be optimal**

Best way to achieve this is by means of automation of slide-making and staining

# Sysmex Staining Solutions



## Semi-Automated Solution

- 5 slides at a time
- Various staining protocols
- Methanol free stain
- 300 slides/kit or 28 days
- Timing, rinsing, drying – controlled.



Kit RAL StainBox MCDh



## Fully Automated Solution

- 10 slides at a time
- Various staining protocols






## Immature Granulocytes – what’s the clinical relevance?



- Manual peripheral blood smear review – includes reporting on presence or absence of “*Left Shift*”

But what does this actually mean to clinicians?

“**Left Shift**” = granulocytic cells in various stages of development.




- The standard Sysmex XN DIFF count includes automated IG count
- The IG count is comprised of metamyelocytes, myelocytes and promyelocytes
- It excludes band cells.
- The presence of IG is **ALWAYS pathological** - with the exception of neonates and in post-partum females ( $\leq 5$ days)



## Clinical utility of the Automated IG count?

- Precursor granulocytes in peripheral blood are pathological
  - Sepsis
  - Haematological disease
  - Bone marrow infiltration
- New cases – suggest use as alert for microscopic review
- Follow-up cases – use automated count for follow-up
  - Sepsis – successful treatment – reduction in IG
  - Bone Marrow recovery post chemotherapy – increase in IG




IG values	Septic Patients	Healthy Donors
Median	0.1	0.01
Minimum	0.0	0.0
Maximum	1.36	0.06
N	40	38

**Automated IG count = enhanced precision and accuracy**  
Useful for detection and monitoring of infections and marrow recovery





# Time course of IG in sepsis patients



## The added value of XN-DIFF

### Features and specifications

added value

**XN-DIFF**

#### **Productivity values:**

- Smear reduction and decreased TAT
- Improved flagging efficiency

#### **Clinical values:**

- High sensitivity detection of WBC abnormalities and efficient flagging
- Support detection of inflammation and/or infection and monitoring response to treatment.

Adding value...

...to the reticulocyte count



### **The XN-RET:**


Providing the complete picture of erythropoiesis using qualitative and quantitative information

- Highly accurate reticulocyte count
- Identification and quantification fragments
- Real-time information on functional iron availability – RET-He



## The Reticulocyte Count – a neglected parameter?

- Most underutilised test in haematology laboratory
- RET# essential for causative diagnosis of anaemia - peripheral destruction or bone marrow production failure?
- Supports early prognosis and optimised treatment serial by monitoring






## RBC sub-population analysis

- Traditional RBC indices provide “average” information of red cell mass produced over past 120 days
- Sub-population analysis provides more dynamic information and aids evaluation of anaemia
- RBC fragments – medical emergency until proven otherwise
  - Automated identification – early warning system


Automated identification of fragments – early warning system for possible for development of microangiopathic haemolysis





## Reticulocyte Haemoglobin Content

- **Real-time** snapshot of haemoglobin concentration in developing RBCs
- Measure of iron availability for erythropoiesis in clinically relevant timeframe



Normal Ret-He



Greatly Reduced Ret-He





## The clinical value of RET-He



A common diagnostic dilemma

Does the patient have true iron deficiency or functional iron deficiency?

Will the patient respond to oral iron supplementation?

# Clinical usefulness of RET-He - cancer patients



American Journal of Clinical Pathology

## Using the Hemoglobin Content of Reticulocytes (RET-He) to Evaluate Anemia in Patients With Cancer



American Society for  
Clinical Pathology

Ellinor I. B. Peerschke, PhD; Melissa S. Pessin, MD, PhD; Peter Maslak, MD

[Disclosures](#) |


Am J Clin Pathol. 2014;142(4):506-512.

- Cancer patients regularly monitored for cytopenias → FBC standard test
- Biochemical iron studies are often included in standard test orders for follow-up visits to speed up management of anaemias
  - RET-he rapidly rules out iron deficiency anaemia
  - Iron studies could be reduced by 80% (209 to 43 in this study)




## Clinical usefulness of RET-He – Renal Patients

- Included in international renal anaemia guidelines to determine when iron supplementation is needed
- Predictor of iron deficiency in haemodialysis patients
- Superior to traditional iron studies



Ferritin: 0.53  
TSAT: 0.56



Ret-He : 0.72

[Buttarello et al, Am J Clin pathol, 2010; 133:949-954](#)



## Clinical usefulness of RET-He



- Reduction in blood transfusion requirements in major orthopaedic surgery – guides when to give EPO to boost pre-operative haemoglobin
- Superior to Ferritin & TSAT in diagnosing iron deficiency – using BM stainable iron as gold standard
- Identification of iron deficiency in children – superior to ferritin and TSAT

# PEDIATRICS®

OFFICIAL JOURNAL OF THE AMERICAN ACADEMY OF PEDIATRICS

### **Diagnosis and Prevention of Iron Deficiency and Iron-Deficiency Anemia in Infants and Young Children (0-3 Years of Age)**


Robert D. Baker, Frank R. Greer and The Committee on Nutrition

*Pediatrics* 2010;126:1040-1050; originally published online Oct 5, 2010;

DOI: 10.1542/peds.2010-2576



## Successful treatment of iron deficiency



RET-He allows for rapid assessment of response to iron supplementation

## The added value of RET: A closer look at erythropoiesis

### Features and specifications

added value

**RET**

#### **Productivity values:**

- Quick screening of erythropoietic status

#### **Clinical values:**

- Quantitative and qualitative data on erythropoietic status of patient to guide anaemia management



Adding value...

...by automating body fluid analysis for rapid diagnosis

added value  
**XN-BF**

**XN-BF:** A dedicated channel for automated cell counting of body fluids with the same precision and accuracy demanded for blood

- Accurate and rapid analysis of body fluids facilitates diagnosis and treatment of several diseases.
- Ability to distinguish between bacterial and non-bacterial infections
- Enhanced accuracy and precision for improved serial monitoring

## CSF analysis – when every second counts

- Question 1: Does the patient have meningitis??

**Answer : ↑ WCC**

Adults > 5-10 cells/ $\mu$ L;

Children > 10-30 cells/ $\mu$ L

- Accurate cell counting at very low levels is vital!
- Question 2: Is it bacterial meningitis??
  - Medical emergency
  - Delayed onset treatment – permanent brain damage or death

**Answer : predominance of polymorphonuclear cells**  
(all other causes of ↑ WCC in CSF are due to mononuclear cells)

- Accurate differential count on body fluid is essential for diagnosis





# Acute Lymphoblastic Leukaemia with meningeal infiltration

## Information from XN-Series

WBC-BF	0.024	10 <sup>9</sup> /L	
RBC-BF	0.000	10 <sup>12</sup> /L	
MN	0.024	10 <sup>9</sup> /L	100.0 %
PMN	0.000	10 <sup>9</sup> /L	0.0 %
TC-BF	0.025	10 <sup>9</sup> /L	

## Research parameters

LY-BF	0.020	10 <sup>9</sup> /L	83.3 %
MO-BF	0.004	10 <sup>9</sup> /L	16.7 %
NE-BF	0.000	10 <sup>9</sup> /L	0.0 %
EO-BF	0.000	10 <sup>9</sup> /L	0.0 %
HF-BF	0.001	10 <sup>9</sup> /L	4.2 /100WBC
RBC-BF2	0.0002	10 <sup>12</sup> /L	



## The added value of XN-BF: the Body Fluids

### Features and specifications

#### **Productivity values:**

- Standardization and reproducibility
- Convenient automatic measurement of a variety of body fluids without need for manual sample preparation
- Rapid TAT

added value

**XN-BF**

#### **Clinical values:**

- Trustworthy results with better precision at clinically relevant cut-off values

Larger laboratories are hooked on XN ...



**Smallest building block –  
100 samples/hour**

The new XN-Series:  
the ideal solution for every s

The standard diagnostic components for each XN module

Analytic Components: XN-CBC and XN-DIFF are standard. All other functi

The XN-1000 and XN-2000 at a glance

**XN-1000**

100 samples/h, 5 rack sampler

200 samples/h

Flexible configurations: from benchtop system to independent workplace Reflex sampler  
rule set Optional wagons for each system for best organization of the workplace reagent

The XN-9000 at a glance

**XN-9000**

<b>Technology</b>	Haematology Automation Analyser	<b>Functional character</b> <ul style="list-style-type: none"><li>• Free configuration of;</li><li>• Discrete rack manage</li><li>• Built-in auto reflex</li><li>• One point of entry/exit</li><li>• Flexible start/stock-yard</li></ul>
<b>Integrated modules</b>	XN analytical, barcode terminal SP for smear making/staining	
<b>Number of modules</b>	max. 9	
<b>Connectable products</b>	tube sorter, ESR, HbA1c	
<b>Throughput</b>	max. 900 samples/h	

XN-Series Information Contact

- ▶ Silent Design
- ▶ Added Value
- ▶ Technology
- ▶ News
- ▶ Contact

Smaller laboratories ...





# The XN-L-Series

---

**XN**

# It truly is XN



- Identical technology
- Flexible model line-up to cater for the needs of any small laboratory



XN-350

Open  
analysis model

XN-450

Closed  
analysis model

XN-550

Sampler  
analysis model

It truly is XN



AL VALUES

PROFESSIONAL SER

PRODUCTIVITY VALUES

- Identical “axis concept”
- Identical concept of “modularity” and “scalability”



# Specifications Portrait



Modes:	WB, PD, L-WBC*, BF*
Aspiration volumes:	25 $\mu$ L in WB and L-WBC mode, 70 $\mu$ L in BF and PD mode Auto dispense function, one aspiration line
Throughput:	60 samples/h, 70 samples/h with Speed-up license*
Parameters:	26 standard diagnostic parameters 15 optional diagnostic parameters
Interface	RS-232C, 2x Ethernet (RJ45), 4x USB
IPU	<b>Integrated</b> , embedded WIN 7 with McAfee Security

\* optional

It truly is XN



- Identical Productivity Value APPs
- Optional flexible speed of the sampler (60 or 70 samples/h)



Rerun & Reflex



Intelligent  
Rack Sampling



Speed-up





L-WBC

It truly is XN



- The same Clinical Value APPs
- Standard: CBC, XN-DIFF
- Optional: RET, XN-BF

# XN-L Series Integrated IPU





## XN-L Series Reagent Concept

- The reagents and the concept for the XN-L Series are the same as for the XN-Series
- The lower throughput of the XN-L has however been taken into consideration
- As a result, smaller reagent pack sizes will be available



## Support for XN and XN-L

- Robust systems – in keeping with the long standing Sysmex reputation
- Installation, end-user training and after sales service support
- Easy transition from XN to XN-L and vice versa
  - Software look and feel the same
  - Reagents can be interchangeable (except RET fluorocell)
  - QC management identical
- Quality management
  - SNCS
  - Sysmex Quality Guidance Manual



XN Concept

